

**Andrew Rader U.S. Army Health Clinic
Fort Myer, VA 22211**



Patient Instructions

PA Form 5

STOOL COLLECTION FOR OVA and PARASITES

Preparation:

1. No special preparations are necessary

Directions:

1. Collect sample onto clean, disposable container.
2. Transfer sufficient sample into each vial to bring the fluid level up to the designated arrow on the label.
3. Write your **Full name and Sponsor's Social Security Number** on each bottle's label.
4. Once the specimens have been placed in the vials, they can be stored at room temperature until you return it to the laboratory
5. Return the vials in the plastic bag provided to the laboratory for testing.